
Product Datasheet

**YWHAZ Antibody (Phospho-Ser58) : FITC
(orb2080905)**

Description

YWHAZ Antibody (Phospho-Ser58) :

Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ELISA, IF, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived between 24-73 amino acids from the Internal region of human 14-3-3 zeta around the phosphorylation site of Ser58.
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Application notes	Application Info: WB: 1:500~1000; IHC: 1:50~100; IF: 1:100~1:500; ELISA: 1:40000. Dilution Info: The optimal working dilution should be determined by the end user.
Clonality	Polyclonal
Uniprot ID	P63104
NCBI	NP_001129171.1
Expiration Date	12 months from date of receipt.