

Product Datasheet

NME1 Antibody - N-terminal region : HRP (orb2081564)

Description	NME1 Antibody - N-terminal region : HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	DOT, IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human NME1
Preservatives	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Form/Appearance	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ANCERTFIAIKPDGVQRGLVGEIIKRFEQKGFRLVGLKFMQASEDLLKEH
Clonality	Polyclonal
MW	17kDa
Uniprot ID	P22392
NCBI	NP_000260
Expiration Date	12 months from date of receipt.