

Product Datasheet

Rab11b Antibody - C-terminal : Biotin (orb2084389)

Description	Rab11b Antibody - C-terminal : Biotin
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	WB
Immunogen	The immunogen for Anti-Rab11b antibody is: synthetic peptide directed towards the C-terminal of Rat RB11B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KNILTEIYRIVSQKQIADRAAHDESPGNNVVDISVPPTTDGQKPNKLQCC
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	O35509
NCBI	NP_116006.1
Expiration Date	12 months from date of receipt.