

---

## Product Datasheet

### Glucose Oxidase antibody (Biotin) (orb20848)

**Description**

Biotin-conjugated IgG fraction of

 oxidase from *Aspergillus ni...*
**Species/Host**

Rabbit

**Conjugation**

Biotin

**Preservatives**

Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Form/Appearance**

Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Storage**

Storage: The lyophilised conjugate is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

**Note**

For research use only

**Clonality**

Polyclonal

**Source**

Glucose oxidase isolated and purified from *Aspergillus niger*. Freund's complete adjuvant is used in the first step of the immunization procedure.

**Uniprot ID**
**P13007**
**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as