

## Product Datasheet

# OR9G9 Antibody - C-terminal region : Biotin (orb2085406)

<b>Description</b>	OR9G9 Antibody - C-terminal region : Biotin
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	WB
<b>Immunogen</b>	The immunogen is a synthetic peptide directed towards the C-terminal region of Human OR9G9
<b>Preservatives</b>	Liquid. Purified antibody supplied in 1x PBS buffer.
<b>Form/Appearance</b>	Liquid. Purified antibody supplied in 1x PBS buffer.
<b>Storage</b>	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Synthetic peptide located within the following region: YIYALPRSSYSFDMDKIVSTFYTEVLPMLNPMIYSLRNKDVKEALKKLLP
<b>Clonality</b>	Polyclonal
<b>MW</b>	33kDa
<b>Uniprot ID</b>	<a href="#">POC7N8</a>
<b>NCBI</b>	<a href="#">NP_001013376</a>
<b>Expiration Date</b>	12 months from date of receipt.