

Product Datasheet

LYPD3 Antibody - C-terminal region : HRP (orb2091146)

Description	LYPD3 Antibody - C-terminal region : HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human LYPD3
Preservatives	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Form/Appearance	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PVRPTSTTKPMPAPTSQTPRQGVHEASRDEEPRLTGGGAAGHQDRSNSGQ
Clonality	Polyclonal
MW	38kDa
Uniprot ID	O95274
NCBI	NP_055215
Expiration Date	12 months from date of receipt.