

Product Datasheet

GGH Antibody - C-terminal region : HRP (orb2093345)

Description	GGH Antibody - C-terminal region : HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	WB
Preservatives	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Form/Appearance	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DGISHAPNAVKTAFYLAEFFVNEARKNNHHFKSESEEEKALIYQFSPIYT
Clonality	Polyclonal
MW	35kDa
Uniprot ID	Q92820
NCBI	NP_003869
Expiration Date	12 months from date of receipt.