

Product Datasheet

SLC25A33 Antibody : FITC (orb2097276)

Description	SLC25A33 Antibody : FITC
Species/Host	Rabbit
Conjugation	FITC
Immunogen	The immunogen is a synthetic peptide directed towards the following sequence SCIAYPHEVIRTRLREEGTKYKSFVQTARLVFREEGYLAFYRGLFAQLIR
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SCIAYPHEVIRTRLREEGTKYKSFVQTARLVFREEGYLAFYRGLFAQLIR
Clonality	Polyclonal
MW	35 kDa
Uniprot ID	Q9BSK2
NCBI	NP_115691
Expiration Date	12 months from date of receipt.