

Product Datasheet

DEFB129 Antibody - C-terminal region : Biotin (orb2100049)

Description	DEFB129 Antibody - C-terminal region : Biotin
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human DEFB129
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TSFFANTNFV IIPNATPMNSATISTMTPGQITYTATSTKSNTKESRDSAT
Clonality	Polyclonal
MW	18kDa
Uniprot ID	Q9H1M3
NCBI	NP_543021
Expiration Date	12 months from date of receipt.