

---

## Product Datasheet

### Alcohol Dehydrogenase antibody (Biotin) (orb21035)

**Description**

Biotin-conjugated IgG fraction of

dehydrogenase from horse li...

**Species/Host**

Rabbit

**Conjugation**

Biotin

**Preservatives**

Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Form/Appearance**

Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Storage**

Storage: The lyophilised conjugate is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

**Note**

For research use only

**Clonality**

Polyclonal

**Source**

Alcohol dehydrogenase isolated and purified from *Leuconostoc mesenteroides*. Freund's complete adjuvant is used in the first step of the immunization procedure.

**Uniprot ID**
**Q5RLY6**
**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO. not for use in diagnostic or