

Product Datasheet

PIK3R5 Antibody - N-terminal region : HRP (orb2105699)

Description	PIK3R5 Antibody - N-terminal region : HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human PIK3R5
Preservatives	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Form/Appearance	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: HRFLTWPVPYCSICQELLTFIDAELKAPGISYQRLVRAEQGLPIRSHRSS
Clonality	Polyclonal
MW	97kDa
Uniprot ID	Q8WYR1
NCBI	NP_055123
Expiration Date	12 months from date of receipt.