

Product Datasheet

RAD9B Antibody - C-terminal region : FITC (orb2108463)

Description	RAD9B Antibody - C-terminal region : FITC
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human RAD9B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ATHAPISYFDYFPGKPLALSIDDMLVEANFILATLADEQSRASSPQSLCL
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q6WBX8
Expiration Date	12 months from date of receipt.