

Product Datasheet

GP1BA Antibody - middle region : FITC (orb2110698)

Description	GP1BA Antibody - middle region : FITC
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human GP1BA
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TPTPKLEKLSLANNLTELPAGLLNGLNLDLLLLQENSLYTIPIKGFSGS
Clonality	Polyclonal
MW	38kDa
Uniprot ID	P07359
NCBI	NP_000164.5
Expiration Date	12 months from date of receipt.