

Product Datasheet

ADCY10 Antibody : Biotin (orb2115682)

Description	ADCY10 Antibody : Biotin
Species/Host	Rabbit
Conjugation	Biotin
Immunogen	The immunogen is a synthetic peptide directed towards the following sequence LMFVDISGFTAMTEKFSSAMYMDRGAEQLVEILNYHISAIVEKVLIFGGD
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LMFVDISGFTAMTEKFSSAMYMDRGAEQLVEILNYHISAIVEKVLIFGGD
Clonality	Polyclonal
MW	187 kDa
Uniprot ID	Q96PN6
NCBI	NP_060887
Expiration Date	12 months from date of receipt.