

Product Datasheet

RAPGEF2 Antibody - C-terminal region : HRP (orb2121695)

Description	RAPGEF2 Antibody - C-terminal region : HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of mouse RAPGEF2
Preservatives	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.
Form/Appearance	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GHTHFDYSGDAASIWASGGHMDQMMFSDHSTKYNRQNSRESLEQAQSR
Clonality	Polyclonal
MW	164kDa
Uniprot ID	Q8CHG7
NCBI	XP_006502317.1
Expiration Date	12 months from date of receipt.