

Product Datasheet

IFIH1 Antibody : HRP (orb2131782)

Description	IFIH1 Antibody : HRP
Species/Host	Rabbit
Conjugation	HRP
Immunogen	The immunogen is a synthetic peptide directed towards the following sequence GSGKTRVAVYIAKDHLDKKKKASEPGKVIVLVNKVLLVEQLFRKEFQPFL
Preservatives	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Form/Appearance	Liquid. Purified antibody is supplied in high phosphate PBS, 100 mm phosphate, 150 mM NaCl, pH 7.6.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GSGKTRVAVYIAKDHLDKKKKASEPGKVIVLVNKVLLVEQLFRKEFQPFL
Clonality	Polyclonal
MW	117 kDa
Uniprot ID	Q9BYX4
NCBI	NP_071451
Expiration Date	12 months from date of receipt.