



Product Datasheet

FGFR1 (phospho-Y654) antibody (orb213927)

Description	Rabbit polyclonal antibody to FGFR1
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y654 of human FGFR1 protein. The exact sequence is proprietary.
Target	FGFR1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q04589 , P16092 , P11362
Entrez	79114 , 14182 , 2260
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of FGFR1 (Phospho-...)



Immunofluorescent analysis of FGFR1 (Pho...)