



Product Datasheet

Kv1.3 (phospho-Y187) antibody (orb214148)


Description	Rabbit polyclonal antibody to KCNA3
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y187 of human Kv1.3 protein. The exact sequence is proprietary.
Target	KCNA3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P22001, P16390, P15384
Entrez	29731, 16491, 3738
Dilution Range	WB: 1-500-1-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Kv1.3 (Phospho-...



Immunohistochemical analysis of Kv1.3 (P...



Immunofluorescent analysis of Kv1.3 (Pho...