


## Product Datasheet

# Nucleophosmin (phospho-T234) antibody (orb214333)

|                            |                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to NPM1                                                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse, Rabbit, Rat                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                          |
| <b>Tested Applications</b> | WB                                                                                                                                                    |
| <b>Immunogen</b>           | KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T234 of human Nucleophosmin protein. The exact sequence is proprietary. |
| <b>Target</b>              | NPM1                                                                                                                                                  |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                             |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                             |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks.<br>For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.            |
| <b>Note</b>                | For research use only                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                                            |
| <b>Source</b>              | Rabbit                                                                                                                                                |
| <b>Uniprot ID</b>          | <b>P13084, P06748, Q61937</b>                                                                                                                         |
| <b>Entrez</b>              | <b>18148, 4869, 25498</b>                                                                                                                             |
| <b>Dilution Range</b>      | WB: 1-500-1-1000                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                       |



Western blot analysis of Nucleophosmin (...)