

## Product Datasheet

# SOSTDC1 ELISA Kit (Rat) [Out of stock] (orb2145133)

<b>Description</b>	SOSTDC1 ELISA Kit (Rat) [Out of stock]
<b>Range</b>	78-5000pg/mL
<b>Tested Applications</b>	ELISA
<b>Storage</b>	Store as indicated in product manual.
<b>Note</b>	For research use only
<b>Application notes</b>	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Uniprot ID</b>	<a href="#">Q642G2</a>
<b>Sensitivity</b>	39 pg/mL
<b>NCBI</b>	<a href="#">NP_714959.1</a>
<b>Expiration Date</b>	Please enquire.