

Product Datasheet

EPHB4 ELISA Kit (Human) (orb2146077)

Description	EPHB4 ELISA Kit (Human)
Range	15.6-1000pg/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, plasma and other biological fluids
Uniprot ID	P54760
Sensitivity	8.04 pg/mL
NCBI	NP_004435.3
Expiration Date	Please enquire.