

Product Datasheet

TNFRSF10C ELISA Kit (Human) [Out of stock] (orb2146563)

Description	TNFRSF10C ELISA Kit (Human) [Out of stock]
Range	0.156-10ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, plasma and other biological fluids
Uniprot ID	O14798
Sensitivity	0.078 ng/mL
NCBI	NP_003832.2
Expiration Date	Please enquire.