

Product Datasheet

NCOA1 ELISA Kit (Human) (orb2147118)

Description	NCOA1 ELISA Kit (Human)
Range	0.156-10ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue homogenates, Cell culture supernates and Other biological fluids.
MW	157 kDa
Uniprot ID	Q15788
Sensitivity	0.056 ng/mL
NCBI	NP_003734.3
Expiration Date	Please enquire.