

Product Datasheet

ANXA5 ELISA Kit (Human) (orb2147388)

Description	ANXA5 ELISA Kit (Human)
Range	31.2-2000pg/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue homogenates, Cell culture supernates, Other biological fluids
MW	36 kDa
Uniprot ID	P08758
Sensitivity	10 pg/mL
NCBI	NP_001145.1
Expiration Date	Please enquire.