

Product Datasheet

IL2 ELISA Kit (Pig) (orb2147610)

Description	IL2 ELISA Kit (Pig)
Range	31.2-2000pg/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue Homogenates , Cell Culture Supernates, Other Biological Fluids
MW	17 kDa
Uniprot ID	P26891
Sensitivity	8.8 pg/mL
NCBI	NP_999026.1
Expiration Date	Please enquire.