

Product Datasheet

PA1 ELISA Kit (Human) (orb2147841)

Description	PA1 ELISA Kit (Human)
Range	0.78-50ng/ml
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Competitive ELISA
Sample Types	Serum, Plasma, Tissue Homogenates, Cell Culture Supernates, Other Biological Fluids
MW	28 kDa
Uniprot ID	Q9BTK6
Sensitivity	0.39 ng/mL
NCBI	NP_078792.1
Expiration Date	Please enquire.