

Product Datasheet

MUC1 ELISA Kit (Human) [Out of stock] (orb2148528)

| | |
|----------------------------|---|
| Description | MUC1 ELISA Kit (Human) [Out of stock] |
| Range | 0.156-10ng/mL |
| Tested Applications | ELISA |
| Storage | Store as indicated in product manual. |
| Note | For research use only |
| Application notes | Assay Info: Assay Methodology: Quantitative Sandwich ELISA |
| Sample Types | Serum, Plasma, Tissue homogenates, Cell culture supernates, Other biological fluids |
| MW | 122 kDa |
| Uniprot ID | P15941 |
| Sensitivity | 0.1 ng/mL |
| NCBI | NP_001018016.1 |
| Expiration Date | Please enquire. |