

## Product Datasheet

### EIF2AK3 ELISA Kit (Human) (orb2149163)

<b>Description</b>	EIF2AK3 ELISA Kit (Human)
<b>Range</b>	0.312-20ng/mL
<b>Tested Applications</b>	ELISA
<b>Storage</b>	Store as indicated in product manual.
<b>Note</b>	For research use only
<b>Application notes</b>	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
<b>Sample Types</b>	Serum, Plasma, Tissue Homogenates , Cell Culture Supernates, Other Biological Fluids
<b>MW</b>	110 kDa
<b>Uniprot ID</b>	<a href="#">Q9NZJ5</a>
<b>Sensitivity</b>	0.15 ng/mL
<b>NCBI</b>	<a href="#">NP_004827.4</a>
<b>Expiration Date</b>	Please enquire.