

Product Datasheet

TNFRSF10B ELISA Kit (Human) [Out of stock] (orb2149429)

Description	TNFRSF10B ELISA Kit (Human) [Out of stock]
Range	0.31-20.0ng/mL
Tested Applications	ELISA
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue Homogenates , Cell Culture Supernates, Other Biological Fluids
MW	48 kDa
Uniprot ID	O14763
Sensitivity	0.16 ng/mL
NCBI	NP_003833.4
Expiration Date	Please enquire.