



Product Datasheet

PAK4 (phospho-S474) antibody (orb214975)


Description	Rabbit polyclonal antibody to PAK4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S474 of human PAK4 protein. The exact sequence is proprietary.
Target	PAK4
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O96013, Q8BTW9
Entrez	10298, 70584
Dilution Range	WB: 1-500-1-1000, IF/ICC: 1-100-1-500
Expiration Date	12 months from date of receipt.



Western blot analysis of PAK4 (Phospho-S474) antibody



Immunohistochemical analysis of PAK4 (Phospho-S474) antibody



Immunofluorescent analysis of PAK4 (Phospho-S474) antibody