

Product Datasheet

IL22RA2 ELISA Kit (Mouse) [Out of stock] (orb2151121)

Description	IL22RA2 ELISA Kit (Mouse) [Out of stock]
Range	0.312-20ng/mL
Tested Applications	ELISA
Storage	12 months shelf life. For unopened kits: All the reagents should be kept at 4C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q80XF5
Sensitivity	0.119ng/mL
Expiration Date	Please enquire.