

## Product Datasheet

### NFKBIA ELISA Kit (Human) (orb2152983)

<b>Description</b>	NFKBIA ELISA Kit (Human)
<b>Range</b>	0.156-10ng/mL
<b>Tested Applications</b>	ELISA
<b>Storage</b>	2°C to 8°C -20°C
<b>Note</b>	For research use only
<b>Sample Types</b>	Tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<a href="#">P25963</a>
<b>Sensitivity</b>	0.065ng/mL
<b>NCBI</b>	<a href="#">NP_065390.1</a>
<b>Expiration Date</b>	Please enquire.