

Product Datasheet

CXCL6 ELISA Kit (Human) (orb2153169)

Description	CXCL6 ELISA Kit (Human)
Range	15.6-1,000pg/mL
Tested Applications	ELISA
Storage	2°C to 8°C -20°C
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P80162
Sensitivity	6.5pg/mL
NCBI	NP_002984.1
Expiration Date	Please enquire.