

Product Datasheet

TNFRSF1A ELISA Kit (Mouse) (orb2153215)

Description	TNFRSF1A ELISA Kit (Mouse)
Range	15.6-1,000pg/mL
Tested Applications	ELISA
Storage	2°C to 8°C -20°C
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P25118
Sensitivity	28pg/mL
NCBI	NP_035739.2
Expiration Date	Please enquire.