

## Product Datasheet

### ANXA5 ELISA Kit (Human) (orb2154233)

<b>Description</b>	ANXA5 ELISA Kit (Human)
<b>Range</b>	31.2-2,000pg/mL
<b>Tested Applications</b>	ELISA
<b>Storage</b>	2°C to 8°C -20°C
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates and other biological fluids
<b>Uniprot ID</b>	<b>P08758</b>
<b>Sensitivity</b>	14.1pg/mL
<b>NCBI</b>	<b>NP_001145.1</b>
<b>Expiration Date</b>	Please enquire.