

## Product Datasheet

# IL-4 ELISA Kit (Guinea Pig ) [Out of stock] (orb2154458)

<b>Description</b>	IL-4 ELISA Kit (Guinea Pig ) [Out of stock]
<b>Range</b>	3.12-200pg/mL
<b>Tested Applications</b>	ELISA
<b>Storage</b>	2°C to 8°C -20°C
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<a href="#">G0WM66</a>
<b>Sensitivity</b>	1.24pg/mL
<b>NCBI</b>	<a href="#">NP_001244192.1</a>
<b>Expiration Date</b>	Please enquire.