

## Product Datasheet

# CXCL6 Chemi-Luminescent ELISA Kit (Rat) (orb2154676)

<b>Description</b>	CXCL6 Chemi-Luminescent ELISA Kit (Rat)
<b>Range</b>	1.37-1,000pg/mL
<b>Tested Applications</b>	CLIA
<b>Storage</b>	2°C to 8°C -20°C
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Uniprot ID</b>	<a href="#">P97885</a>
<b>Sensitivity</b>	0.49pg/mL
<b>NCBI</b>	<a href="#">NP_071550.1</a>
<b>Expiration Date</b>	Please enquire.