

Product Datasheet

IGFBP4 Chemi-Luminescent ELISA Kit (Human) [Out of stock] (orb2155425)

Description	IGFBP4 Chemi-Luminescent ELISA Kit (Human) [Out of stock]
Range	0.55-400ng/mL
Tested Applications	CLIA
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P22692
Sensitivity	0.55 ng/mL
Expiration Date	Please enquire.