

Product Datasheet

HNF1A Chemi-Luminescent ELISA Kit (Mouse) [Out of stock] (orb2155467)

Description	HNF1A Chemi-Luminescent ELISA Kit (Mouse) [Out of stock]
Range	13.7-10,000pg/mL
Tested Applications	CLIA
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Tissue homogenates and other biological fluids
Uniprot ID	P22361
Sensitivity	0.156 ng/mL
NCBI	NP_033353
Expiration Date	Please enquire.