

Product Datasheet

CXCL6 Chemi-Luminescent ELISA Kit (Human) [Out of stock] (orb2155519)

Description	CXCL6 Chemi-Luminescent ELISA Kit (Human) [Out of stock]
Range	2.74-2,000pg/mL
Tested Applications	CLIA
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, plasma and other biological fluids
Uniprot ID	P80162
Sensitivity	31.25 pg/mL
NCBI	NP_002984
Expiration Date	Please enquire.