

## Product Datasheet

### LTA4H ELISA Kit (Human) (orb2156800)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | LTA4H ELISA Kit (Human)   |
| <b>Range</b>               | 0.312-20ng/mL   |
| <b>Tested Applications</b> | ELISA   |
| <b>Storage</b>             | 2°C to 8°C -20°C  |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Assay Info: Assay Methodology:<br>Quantitative Sandwich Immunoassay |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates and other biological fluids       |
| <b>Uniprot ID</b>          | <b>P09960</b>   |
| <b>Sensitivity</b>         | 0.128 ng/mL   |
| <b>NCBI</b>                | <b>NP_000886.1</b>  |
| <b>Expiration Date</b>     | Please enquire.   |