

Product Datasheet

IL22RA2 ELISA Kit (Rat) (orb2157481)

Description	IL22RA2 ELISA Kit (Rat)
Range	0.312-20ng/mL
Tested Applications	ELISA
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich Immunoassay
Sample Types	Serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q7TNI4
Sensitivity	0.105 ng/mL
NCBI	NP_001003404
Expiration Date	Please enquire.