

Product Datasheet

CD9 ELISA Kit (Human) [Out of stock] (orb2157670)

Description	CD9 ELISA Kit (Human) [Out of stock]
Range	31.2-2,000pg/mL
Tested Applications	ELISA
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich Immunoassay
Sample Types	Serum, plasma, tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P21926
Sensitivity	12.5 pg/mL
NCBI	NP_001760
Expiration Date	Please enquire.