

Product Datasheet

EXOSC2 ELISA Kit (Human) [Out of stock] (orb2157871)

Description	EXOSC2 ELISA Kit (Human) [Out of stock]
Range	0.156-10ng/mL
Tested Applications	ELISA
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich Immunoassay
Sample Types	Tissue homogenates, cell lysates and other biological fluids
Uniprot ID	Q13868
Sensitivity	0.057 ng/mL
NCBI	NP_001269637
Expiration Date	Please enquire.