

Product Datasheet

SOSTDC1 ELISA Kit (Rat) (orb2158318)

Description	SOSTDC1 ELISA Kit (Rat)
Range	23.5 - 1,500 pg/mL
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich Immunoassay
Sample Types	Serum, Plasma, Tissue Homogenates
Uniprot ID	Q642G2
Sensitivity	5.86 pg/mL
NCBI	NP_714959
Expiration Date	Please enquire.