

Product Datasheet

TNFRSF10D ELISA Kit (Human) [Out of stock] (orb2158412)

Description	TNFRSF10D ELISA Kit (Human) [Out of stock]
Range	0.312 - 20 ng/mL
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich Immunoassay
Sample Types	Serum, Plasma, Tissue Homogenates
Uniprot ID	Q9UBN6
Sensitivity	0.078 ng/mL
NCBI	NP_003831
Expiration Date	Please enquire.