


Product Datasheet

Human Kininogen-1 / KNG1 / HMW Kininogen ELISA Kit (orb215856)

Description	Human Kininogen-1 / KNG1 / HMW Kininogen ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q19-S644
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:500, concentration in serum and plasma is around 600 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	71957 MW
Uniprot ID	P01042
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human Kininogen-1/KNG1 ELISA Kit standard...