

Product Datasheet

IGF2BP2 ELISA Kit (Human) [Out of stock] (orb2158587)

Description	IGF2BP2 ELISA Kit (Human) [Out of stock]
Range	125 - 8,000 pg/mL
Storage	Store as indicated in product manual.
Note	For research use only
Application notes	Assay Info: Assay Methodology: Quantitative Sandwich ELISA
Sample Types	Serum, Plasma, Tissue Homogenates, Cell Lysates
Uniprot ID	Q9Y6M1
Sensitivity	31.25 pg/mL
NCBI	NP_001007226.1
Expiration Date	Please enquire.