


Product Datasheet

Human TFPI2 ELISA Kit (orb215864)

Description	Human TFPI2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D23-F235
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 2 ng/ml
Sample Types	cell culture supernatants and serum.
MW	26934 MW
Uniprot ID	P48307
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human TFPI2
PicoKine
ELISA Kit
standard ...