


Product Datasheet

Human ADAMTS4 ELISA Kit (orb215872)

Description	Human ADAMTS4 ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: F213-C685
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 625 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	90197 MW
Uniprot ID	075173
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human ADAMTS4 PicoKine ELISA Kit standar...